

For more information on **Directions for accessing PICC line using Intermate® in the home setting**, please contact Cathy Hopkins CNS in Our Lady's Children's Hospital, Crumlin, Dublin 12.

Contact details office hours:

- Phone 4096100 Bleep 543/523
- Admin office 4096096/6338/6893

Contact details out of hours:

- 01-4096100 main switch and ask to be put through to the Infectious Disease consultant on call

Additional instructions

PICC Size _____

Date of Insertion _____

Site _____

Remember: Ask your nurse if you are unsure about anything about your child's care.

Developed by: Catherine Hopkins OPAT CNSp, Dr. Patrick Gavin, Amanda Walsh, OPAT Team OLCHC

Date issued: Nov 2015

Date of review: Nov 2018

©2015, Our Lady's Children's Hospital Crumlin, Dublin 12. All rights reserved. No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means without the prior written permission of the copyright holder. Every effort has been made to ensure that the information provided is accurate and in accord with standards accepted at the time of printing.

Our Lady's Children's Hospital, Crumlin

Our Lady's Children's Hospital, Crumlin, Dublin 12

...where children's health comes first

Information for parents/carers of a child on Outpatient Antibiotic Therapy (OPAT)

Directions for accessing PICC line using Intermate® device in the home setting

(1Fr Cook/Vygon. 2Fr Cook/Vygon)

1. Preparing Area:

- Wash hands as per leaflet provided
- Clean tray or table with wipes provided by Baxter
- Open pack on to a cleaned area
- Check details on Intermate® match prescription
 - Correct child's name
 - Correct antibiotic name and dose
 - Correct date and time
- Place Intermate® on prepared area
- Wash hands again as per leaflet provided
- Open syringes by taking off white caps. Remove any air. Replace white caps.
- Place on prepared area.
- Open alcohol swabs on the plastic glove (as they dry quickly)

2. Priming Intermate®

- Remove blue winged cap at the end of the tubing
- Open white clamp
- Let 3-4 drops of fluid out. Line is now primed of air
- Replace blue winged cap
- Close white clamp

3. Cleaning needle free device, non-touch technique:

- Wash hands as per leaflet
- Open alcowipe, hold in non-dominant hand, use to pick up PICC
- Open second alcowipe and firmly scrub bionector for 30 seconds
- Let dry for 20 seconds
- Attach Saline flush NACL 0.9% syringe, pushing hard and turning to the right
- Open PICC clamp
- Flush in 3mls Saline flush NACL 0.9% using the push pause technique (to create turbulence)
- Close PICC clamp.
- Remove blue winged cap
- Attach antibiotic (push in hard and to the right)
- Open all clamps
- Close over pack and keep clean for use after antibiotic infused.

4. Ensure proper infusion

- Check tubing is not kinked and clamps are open
- "Balloon" will deflate over ____ mins approx.
- Infusion complete when "balloon" completely deflated and 8 indicator bumps are visible

5. After Antibiotic infused

- Clamp PICC
- Unscrew antibiotic (Intermate device)
- Wash hands as per leaflet
- Open pack
- Open alcowipe, hold alcowipe in non-dominant hand, use to pick up PICC
- Pick up another alcowipe and firmly scrub needle free device for 30 seconds.
- Let dry for 20 seconds
- Attach Saline flush NACL 0.9% pushing in hard and turning to the right
- Open PICC clamp
- Flush 3mls quickly using the push pause technique
- Close PICC clamp
- Attach heparin sodium syringe (10 units/ml) anticoagulant to stop clotting.
- Open PICC clamp
- Push in _____ml, closing clamp while pushing in another _____ml (positive pressure technique).
- Close PICC clamp
- Dispose Intermate® and syringe in the treatment waste bin provided