

Our Lady's Children's Hospital Community Care Request for Oxygen Saturation Monitor & Supplies

Patient Details	Diagnosis	
	Weight	
Parent / Guardian:	Oxygen Saturation Parameters	
Contact Number:	Upper Lower	
	Heart Data Davameters	
	Heart Rate Parameters	
	Upper Lower	
Medical Card No :	Saturation Monitoring Continuous	
	Spot check	
Disch anging Mand	Tall	
Discharging Ward : Clinical Nurse Specialist :	Tel : Tel :	
Discharge Coordinator:	Tel:	
F	lease Tick as appropriate	
Saturation Monitor 1 unit		
2 units		
If Continuous monitoring is required 2 units should be provided with one being hand held.		
Saturation Probe Sited Finger		
Foot		
Toe [
Approximate Monthly Probe Use	per month	
Prescribing Hospital :	IMO	
Prescribers Signature : Prescribers Name :	IMC:	
FICSUIDEIS INGILIE.		

Issued: February 2015, Review February 2016 Duplicate Copy to be filed in Health Care Records



Oxygen Saturation Monitors provided should meet the following specifications. A list of suppliers is attached, OLCHC are unable to provide any costing for equipment or related consumables.

Specification:

SPO2: 1%-100% [+/- 1%]

Heart Rate: 20 to 250 beats per minute

Audible and Visual alarms for: High/Low Saturation

High/Low Pulse Rate

Sensor Off

Low Battery

Features:

User configurable power up default settings

Adjustable Audible Alarm Volumes

Protective boot cover (If handheld)

Power:

Interchangeable batteries, e.g. AA (Mains operation with minimum of 4hrs Battery Backup @ 100% Charge)

Display:

%SPO2, Heart Rate

Visual indicators for: Pulse search

Audible alarms silenced or off

Battery capacity

Suppliers

Companies will advise on the appropriate model and probes according to the individual monitoring needs of each patient as listed above.

Monitor Name	Supplier	Contact
Nellcor	Health Care 21	01 885 5551
Massimo	Norso Medical	+44 2890343927
		www.norsomedical.com
BCI SPO2	Smiths Medical	B Maguire
		01 294 1133

Issued: February 2015, Review February 2016 Duplicate Copy to be filed in Health Care Records